

In [3]:

```
import numpy as np
import pandas as pd
import os
import nibabel as nib
import pydicom as dicom
import matplotlib.pyplot as plt
```

## Convert nii file to txt file and plot ¶

In [4]:

```
HOME = 'D:\\results0622\\Results\\'
subdirs = ['ROISignals_FunImgWCF']
```

In [5]:

```
for subdir in subdirs:
    home = HOME + subdir
    for fname in os.listdir(home):
        if fname.split('.')[0] == 'nii':
            fp = home + '\\\\' + fname
            f = nib.load(fp)
            f_data = f.get_fdata()
            dim_0, dim_1, dim_2 = f.header.get_data_shape()
            N = dim_0 * dim_1 * dim_2
            np.savetxt(home + '\\\\' + fname.split('.')[0] + '.txt',
                       f_data.reshape(N,))
```

In [16]:

```
# a = np.array([
#
#         [[0,1,2],
#          [3,4,5]],
#
#         [[6,7,8],
#          [9,10,11,]],
#
#         [[12,13,14],
#          [15,16,17]],
#
#         [[18,19,20],
#          [21,22,23]]
#
#         ])
# print(a.shape)
# print(a.reshape(24,))
```


```
(4, 2, 3)
[ 0  1  2  3  4  5  6  7  8  9 10 11 12 13 14 15 16 17 18 19 20 21 22 23]
```

In [16]:

```
subjects = ['Sub_001', 'Sub_002', 'Sub_003', 'Sub_004', 'Sub_005', 'Sub_006', 'Sub_007', 'Sub_008', 'Sub_009', 'Sub_010',]  
p = r'D:\results0622\Results\ROISignals_FunImgWCF'
```


In [17]:

```
for i in subjects:  
    matrix = np.loadtxt(p + '\\ROICorrelation_FisherZ_'+ i +'.txt')  
    plt.imshow(matrix)  
    plt.title('ROICorrelation_FisherZ_' + i)  
    plt.savefig(p + '\\ROICorrelation_FisherZ_' + i +'.png')
```




In [18]:

```
for i in subjects:  
    matrix = np.loadtxt(p + '\\ROICorrelation_'+ i +'.txt')  
    plt.imshow(matrix)  
    plt.title('ROICorrelation_' + i)  
    plt.savefig(p + '\\ROICorrelation_' + i +'.png')
```



In [19]:

```
for i in subjects:
    matrix = np.loadtxt(p + '\\ROISignals_'+ i +'.txt')
    plt.imshow(matrix)
    plt.title('ROISignals_' + i)
    plt.savefig(p + '\\ROISignals_' + i +'.png')
```



In [89]:

```

for i in subjects:
    p = r'D:\result0620\Results\FunImgARFCB_ROISignals\ROISignals_' + i + '.txt'
    s = np.loadtxt(p)
    print('ROISignals_' + i + ': ', s.shape)
for i in subjects:
    p = r'D:\result0620\Results\FunImgARFCB_ROISignals\ROICorrelation_' + i + '.txt'
    s = np.loadtxt(p)
    print('ROICorrelation_' + i + ': ', s.shape)
for i in subjects:
    p = r'D:\result0620\Results\FunImgARFCB_ROISignals\ROICorrelation_FisherZ_' + i +
    '.txt'
    s = np.loadtxt(p)
    print('ROICorrelation_FisherZ_' + i + ': ', s.shape)

```

```

ROISignals_Sub_001: (83, 116)
ROISignals_Sub_002: (107, 116)
ROISignals_Sub_003: (124, 116)
ROISignals_Sub_004: (121, 116)
ROISignals_Sub_005: (89, 116)
ROISignals_Sub_006: (86, 116)
ROISignals_Sub_007: (130, 116)
ROISignals_Sub_008: (116, 116)
ROISignals_Sub_009: (116, 116)
ROISignals_Sub_010: (91, 116)
ROICorrelation_Sub_001: (116, 116)
ROICorrelation_Sub_002: (116, 116)
ROICorrelation_Sub_003: (116, 116)
ROICorrelation_Sub_004: (116, 116)
ROICorrelation_Sub_005: (116, 116)
ROICorrelation_Sub_006: (116, 116)
ROICorrelation_Sub_007: (116, 116)
ROICorrelation_Sub_008: (116, 116)
ROICorrelation_Sub_009: (116, 116)
ROICorrelation_Sub_010: (116, 116)
ROICorrelation_FisherZ_Sub_001: (116, 116)
ROICorrelation_FisherZ_Sub_002: (116, 116)
ROICorrelation_FisherZ_Sub_003: (116, 116)
ROICorrelation_FisherZ_Sub_004: (116, 116)
ROICorrelation_FisherZ_Sub_005: (116, 116)
ROICorrelation_FisherZ_Sub_006: (116, 116)
ROICorrelation_FisherZ_Sub_007: (116, 116)
ROICorrelation_FisherZ_Sub_008: (116, 116)
ROICorrelation_FisherZ_Sub_009: (116, 116)
ROICorrelation_FisherZ_Sub_010: (116, 116)

```

## Version check

In [10]:

```
dicom.__version__
```

Out[10]:

```
'2.0.0'
```

In [11]:

```
nib.__version__
```

Out[11]:

```
'3.1.0'
```

## NII file

In [20]:

```
path = r'D:\results0622\Results\fALFF_FunImgWC'
img = nib.load(path + '\\zfALFFMap_Sub_001.nii')
img_data = img.get_fdata()
print(img.header.get_data_shape())
img_data[6, 29, 21]
```

```
(61, 73, 61)
```

Out[20]:

```
0.6847445964813232
```

In [22]:


```
def show_slices(slices):
    fig, axes = plt.subplots(1, len(slices))
    for i, slice in enumerate(slices):
        axes[i].imshow(slice.T, cmap="gray", origin="lower")

slice_0 = img_data[6, :, :]
slice_1 = img_data[:, 29, :]
slice_2 = img_data[:, :, 21]
show_slices([slice_0, slice_1, slice_2])
plt.suptitle("Center slices for EPI image")
```

Out[22]:

```
Text(0.5, 0.98, 'Center slices for EPI image')
```

Center slices for EPI image



In [ ]:

# DCM file

In [13]:

```
img_hcy = dicom.dcmread(path + '\\FunRaw\\Sub_001\\000001.dcm')  
print(img_hcy)
```

```

Dataset.file_meta -----
(0002, 0000) File Meta Information Group Length  UL: 210
(0002, 0001) File Meta Information Version       OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID        UI: MR Image Storage
(0002, 0003) Media Storage SOP Instance UID     UI: 1.3.6.1.4.1.9590.100.
1.2.374667898811965955833216997681997868070
(0002, 0010) Transfer Syntax UID                UI: Explicit VR Little En
dian
(0002, 0012) Implementation Class UID          UI: 1.3.6.1.4.1.9590.100.
1.3.100.9.4
(0002, 0013) Implementation Version Name       SH: 'MATLAB IPT 9.4'
-----
(0008, 0005) Specific Character Set             CS: 'ISO_IR 100'
(0008, 0008) Image Type                        CS: ['ORIGINAL', 'PRIMAR
Y', 'M_FFE', 'M', 'FFE']
(0008, 0013) Instance Creation Time            TM: '160318'
(0008, 0014) Instance Creator UID              UI: 2.16.124.113543.6006.
99.3051493601938435224
(0008, 0016) SOP Class UID                     UI: MR Image Storage
(0008, 0018) SOP Instance UID                  UI: 1.3.6.1.4.1.9590.100.
1.2.374667898811965955833216997681997868070
(0008, 0020) Study Date                       DA: '20121114'
(0008, 0021) Series Date                      DA: '20121114'
(0008, 0022) Acquisition Date                  DA: '20121114'
(0008, 0023) Content Date                     DA: '20121114'
(0008, 0030) Study Time                       TM: '151828.000000'
(0008, 0031) Series Time                      TM: '153744.79000'
(0008, 0032) Acquisition Time                  TM: '153744.79'
(0008, 0033) Content Time                     TM: '153744.79'
(0008, 0060) Modality                         CS: 'MR'
(0008, 0070) Manufacturer                     LO: 'Philips Medical Syst
ems'
(0008, 0080) Institution Name                  LO: 'OHSU MRI3'
(0008, 0090) Referring Physician's Name       PN: ''
(0008, 1030) Study Description                 LO: 'MRI BRAIN W/O CONTRA
ST'
(0008, 1032) Procedure Code Sequence 1 item(s) ----
  (0008, 010b) Context Group Extension Flag    CS: 'N'
  -----
  (0008, 103e) Series Description              LO: 'Resting State fMRI'
  (0008, 1060) Name of Physician(s) Reading Study PN: ''
  (0008, 1070) Operators' Name                 PN: ''
  (0008, 1080) Admitting Diagnoses Description LO: ''
  (0008, 1090) Manufacturer's Model Name       LO: 'Intera'
  (0008, 1110) Referenced Study Sequence 1 item(s) ----
    (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.06282883507200984267
    (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.6893832446125167224
    -----
  (0008, 1111) Referenced Performed Procedure Step Sequence 1 item(s) ---
  -
    (0008, 0013) Instance Creation Time        TM: '160323'
    (0008, 0014) Instance Creator UID          UI: 2.16.124.113543.60
06.99.3051493601938435224
    (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.06279230344779385614
    (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.03779990973132413033
    (0020, 0013) Instance Number               IS: "0"
    (2005, 1406) Private tag data              UN: b'\x00\x00'

```



```
-----
(0008, 1120) Referenced Patient Sequence    1 item(s) ----
  (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.06278070223545294329
  (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.4513944610792084839
-----
(0008, 1140) Referenced Image Sequence     3 item(s) ----
  (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.4332567756067262732
  (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.03069736621497723066
-----
  (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.4332567756067262732
  (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.03101131755547137486
-----
  (0008, 1150) Referenced SOP Class UID      UI: 2.16.124.113543.60
06.99.4332567756067262732
  (0008, 1155) Referenced SOP Instance UID   UI: 2.16.124.113543.60
06.99.03027800815908766954
-----
(0010, 0010) Patient's Name                 PN: '00000001'
(0010, 0020) Patient ID                     LO: '00000001'
(0010, 0030) Patient's Birth Date           DA: ''
(0010, 0040) Patient's Sex                  CS: 'F'
(0010, 1010) Patient's Age                   AS: '075Y'
(0010, 1030) Patient's Weight                DS: "80.0"
(0010, 21c0) Pregnancy Status                US: 4
(0018, 0015) Body Part Examined              CS: 'BRAIN'
(0018, 0020) Scanning Sequence              CS: 'GR'
(0018, 0021) Sequence Variant               CS: 'SK'
(0018, 0022) Scan Options                   CS: 'FS'
(0018, 0023) MR Acquisition Type             CS: '2D'
(0018, 0024) Sequence Name                  SH: ''
(0018, 0050) Slice Thickness                 DS: "3.313"
(0018, 0080) Repetition Time                 DS: "3000.99829101562"
(0018, 0081) Echo Time                      DS: "30.001"
(0018, 0083) Number of Averages              DS: "1.0"
(0018, 0084) Imaging Frequency              DS: "127.778973"
(0018, 0085) Imaged Nucleus                 SH: '1H'
(0018, 0086) Echo Number(s)                 IS: "1"
(0018, 0087) Magnetic Field Strength         DS: "3.0"
(0018, 0088) Spacing Between Slices         DS: "3.313"
(0018, 0089) Number of Phase Encoding Steps IS: "63"
(0018, 0091) Echo Train Length              IS: "59"
(0018, 0093) Percent Sampling                DS: "100.0"
(0018, 0094) Percent Phase Field of View    DS: "93.75"
(0018, 0095) Pixel Bandwidth                 DS: "2110.0"
(0018, 1000) Device Serial Number            LO: '4sEr9bzYud6E'
(0018, 1020) Software Versions               LO: ['3.2.1', '3.2.1.1']
(0018, 1030) Protocol Name                   LO: 'RestingStatefMRI'
(0018, 1060) Trigger Time                   DS: "0.0"
(0018, 1081) Low R-R Value                   IS: "0"
(0018, 1082) High R-R Value                 IS: "0"
(0018, 1083) Intervals Acquired             IS: "0"
(0018, 1084) Intervals Rejected             IS: "0"
(0018, 1088) Heart Rate                     IS: "0"
(0018, 1100) Reconstruction Diameter         DS: "212.0"
(0018, 1250) Receive Coil Name              SH: 'SENSE-Head-8'
```

```

(0018, 1251) Transmit Coil Name          SH: 'B'
(0018, 1310) Acquisition Matrix          US: [64, 0, 0, 63]
(0018, 1312) In-plane Phase Encoding Direction CS: 'COL'
(0018, 1314) Flip Angle                  DS: "80.0"
(0018, 1316) SAR                         DS: "0.06301491707563"
(0018, 1318) dB/dt                       DS: "45.3371336787068"
(0018, 5100) Patient Position            CS: 'HFS'
(0018, 9073) Acquisition Duration        FD: 423.1409912109375
(0018, 9087) Diffusion b-value           FD: 0.0
(0018, 9089) Diffusion Gradient Orientation FD: [0.0, 0.0, 0.0]
(0020, 000d) Study Instance UID          UI: 2.16.124.113543.6006.
99.6893832446125167224
(0020, 000e) Series Instance UID        UI: 2.16.124.113543.6006.
99.03075240958158977453
(0020, 0010) Study ID                    SH: '16904854'
(0020, 0011) Series Number               IS: "501"
(0020, 0012) Acquisition Number          IS: "5"
(0020, 0013) Instance Number             IS: "1"
(0020, 0032) Image Position (Patient)    DS: [-104.74581944942, -1
02.64433813095, -47.136009216308]
(0020, 0037) Image Orientation (Patient) DS: [1, 0, 0, 0, 1, 0]
(0020, 0052) Frame of Reference UID      UI: 2.16.124.113543.6006.
99.02322024500391392010
(0020, 0100) Temporal Position Identifier IS: "1"
(0020, 0105) Number of Temporal Positions IS: "140"
(0020, 1040) Position Reference Indicator LO: ''
(0020, 1041) Slice Location              DS: "0.0"
(0028, 0002) Samples per Pixel           US: 1
(0028, 0004) Photometric Interpretation  CS: 'MONOCHROME2'
(0028, 0010) Rows                        US: 64
(0028, 0011) Columns                     US: 64
(0028, 0030) Pixel Spacing               DS: [3.3125, 3.3125]
(0028, 0100) Bits Allocated              US: 16
(0028, 0101) Bits Stored                  US: 16
(0028, 0102) High Bit                     US: 15
(0028, 0103) Pixel Representation        US: 0
(0028, 0106) Smallest Image Pixel Value  US: 0
(0028, 0107) Largest Image Pixel Value   US: 1638
(0028, 1050) Window Center                DS: "865.25"
(0028, 1051) Window Width                 DS: "1504.09"
(0028, 2110) Lossy Image Compression      CS: '00'
(0032, 000a) Study Status ID              CS: 'READ'
(0032, 000c) Study Priority ID            CS: 'LOW'
(0040, 0241) Performed Station AE Title   AE: 'MRI3'
(0040, 0245) Performed Procedure Step Start Time TM: '151828'
(0040, 0251) Performed Procedure Step End Time TM: '151828'
(0040, 0260) Performed Protocol Code Sequence 1 item(s) ----
(0008, 010b) Context Group Extension Flag CS: 'N'
-----
(0040, 0275) Request Attributes Sequence 1 item(s) ----
-----
(0040, 0321) Film Consumption Sequence 0 item(s) ----
(0040, 9096) Real World Value Mapping Sequence 1 item(s) ----
(0040, 9224) Real World Value Intercept    FD: 0.0
(0040, 9225) Real World Value Slope        FD: 1.2366300366300367
-----
(2050, 0020) Presentation LUT Shape        CS: 'IDENTITY'
(7fe0, 0010) Pixel Data                    OW: Array of 8192 element
S

```

In [ ]:

In [ ]:

In [ ]: